



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Choline Oxidase from *Alcaligenes sp.*

Product Number **C 5896**
Storage Temperature -0 °C

Product Description

Enzyme Commission (EC) Number: 1.1.3.17
CAS Number: 9028-67-5
Molecular Weight: 66 kDa (SDS-PAGE) and 72 kDa
(gel filtration)¹
pI: 4.1 ± 0.1.

The enzyme can be used for the enzymatic determination of phospholipids by coupling with phospholipase D and for cholinesterase activity assays.^{2,3}

This enzyme has several inhibitors:
p-chloromercuribenzoate, Cu, Co, Hg, and Ag ions.
The optimal pH is 8.0 - 8.5. The optimal temperature is 40 - 45 °C.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Storage/Stability

Solution stability was measured as a plot of activity vs. time for an enzyme concentration of 1.0 mg/ml in 0.1 M potassium phosphate buffer, pH 7.5, at 37 °C. Approximately 75% of the enzymatic activity remained after 16 hours. In the presence of 10 mM EDTA or 0.5 mg/ml BSA, there is approximately 90% enzymatic activity remaining. Addition of both EDTA and BSA resulted in nearly 100% enzymatic activity remaining after 16 hours.

References

1. Ohta-Fukuyama, M., Identification and Properties of the Prosthetic Group of Choline Oxidase from *Alcaligenes Sp.* J. Biochem., **88**, 197-203 (1980).
2. Takayama, M. et al., Clin. Chim. Acta, 79, 93 (1977).
3. Okabe, H. et al., Clin. Chim. Acta, 80, 87 (1977).

MWM/RXR 9/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.